

## Nitro HD Rollator

HCPCS:




FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Attractive, Euro-style design
- Brake cable inside frame for added safety (A)
- Handle height easily adjusts with unique push button
- Unique suspension system for support and comfort (B)
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius (C)
- Large 10" front casters allow optimal steering and rolling comfort (C)
- Easily folds with one hand to ultra-compact size for storage (D)
- Under seat locking mechanism holds the Nitro closed when frame is folded (E)
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable (F)
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight, aluminum frame

### PRODUCTS

	Description
10266HD-BK	HD Nitro Rollator, Black, 1/ea

### SPECIFICATIONS

	Description
Handle Height	36" - 41"
Seat Dimensions	23" (W) x 10" (D)
Seat Height	23.5"
Unit Dimensions	28" (L) x 28.25" (W)
Warranty	Frame: Limited Lifetime   Brake Handle and Housing: 5 Years
Weight	23 lbs
Weight Capacity	450 lbs

### ACCESSORIES

	Description
10266-BRAKE ENH	Nitro Slow Down Brake, ENH, 1PR, 25/CS, EA

### PARTS

	Description
1026646	CrossFrameBolt,Nitro/NitroHemiHD/NitroH, EA
10266HD-33	Cross Frame, Nitro HD, 1/ea, EA
10266HD01-R	Hand Brake Assmby, Right, Nitro HD, 1ea, EA
1026617-AN	FrontWheelW/Bearings,NitroW/NewFork,1/e, EA
1026611	ToteBagMountSetW/Screws for Nitro, 1/pr, PR
1026641	WheelBrakeLeftW/Hrdwr,Plastic,Nitro,1ea, EA

PARTS	Description
1026625-AN	Wheel&ForkAssbly,Left,Alum,NitroNew,1/ea, EA
1026624-AN	Wheel&ForkAssbly,Rght,Alum,NitroNew,1/ea, EA
10266HD01-L	Hand Brake Assmbly, Left, Nitro HD, 1/ea, EA
1026642	WheelBrakeRightW/Hrdwr,Plastic,Nitro,1e, EA
1026644	Folding Lock Mechanism, Red, Nitro,1/ea, EA
1026618-AN	FrontWheelBearings,NitroW/NewFork, 1/pr, PR
1026614-AN	Hub Cap for Aluminum NitroAlum.Fork,1/ea, EA
1026622-AN	Fork, Left, Aluminum, Nitro New, 1/ea, EA
1026643	Cross Brace Screw & Insert, Nitro, 1/ea, EA
1026621-AN	Fork, Right, Aluminum, Nitro New, 1/ea, EA
1026645	Front Axle & Washer,NitroW/NewFork, 1/e, EA
892	Cane Holder, Nitro Rollator, 1/ea, EA
1026637	Step Tube, Right, Nitro, 1/ea, EA
10266HD31	Shock Absorber, Nitro HD, 1/ea, EA