

# Durable medical equipment **product guide**



Innovative products and support




# Elevate care with medical equipment that goes the distance

Find lots of new products inside including folding options for bath safety items and additional arm and seat options for wheelchairs

## What's inside:

Bath safety.....	3
Wheelchairs.....	11
Transport chairs.....	18
Rollators.....	20
Walkers.....	26
Knee walkers.....	30
Homecare beds & accessories.....	31
Canes.....	34



# Bath safety

## 3-in-1 commodes

- Use as a bedside commode or over toilet
- Folding commode reduces storage space by 35%

Item No.	Description	Wt. capacity	Seat height	Seat depth	Inside arm width	Overall width	Pkg.
MDS89664	Grey, fixed, standard seating	350 lbs. (158.8 kg)	15"-21"	15"	19"	24"	4 cs
G1-200X1	White, folding, standard seating	350 lbs. (158.8 kg)	16"-22"	14.5"	19"	21"	1 ea
G1-200X4	White, folding, standard seating	350 lbs. (158.8 kg)	16"-22"	14.5"	19"	21"	4 cs
G1-200LX1	White, folding, elongated seating	350 lbs. (158.8 kg)	16"-22"	16.5"	19"	21"	1 ea
G1-200LX4	White, folding, elongated seating	350 lbs. (158.8 kg)	16"-22"	16.5"	19"	21"	4 cs

## Drop arm commodes

- Armrests conveniently swing away for easy patient transfers; padded option available for increased comfort

Item No.	Description	Wt. capacity	Seat height	Seat depth	Inside arm width	Overall width	Pkg.
G1-301DX1	White, drop arm, elongated	300 lbs. (136.1 kg)	19"-23"	16.5"	20"	23"	1 ea
G1-301DX2	White, drop arm, elongated	300 lbs. (136.1 kg)	19"-23"	16.5"	20"	23"	2 ea
G1-301DPX1	White, padded, drop arm, elongated	300 lbs. (136.1 kg)	19"-23"	16.5"	20"	23"	1 ea
G1-303DWX1	White drop arm elongated, 4 locking casters	300 lbs. (136.1 kg)	17"-22.5"	17.5"	16.5"	23"	1 ea

Drop arm with 300 lbs. weight capacity



G1-200X1



G1-200LX1



G1-301DX1



G1-303DWX1

## Economy commodes

- Basic commode without lid or splash guard; available with standard or elongated seat

Item No.	Description	Wt. capacity	Seat height	Seat depth	Inside arm width	Overall width	Pkg.
G1-202EX4	Standard steel folding commode	350 lbs. (158.8 kg)	16"-22"	14.5"	19"	21"	4 cs
G1-202ELX4	Elongated steel folding commode	350 lbs. (158.8 kg)	16"-22"	16.5"	19"	21"	4 cs

## Bariatric commodes

- Fixed and folding commodes designed with sturdy steel frames for high weight capacity needs

Item No.	Description	Wt. capacity	Seat height	Seat depth	Inside arm width	Overall width	Pkg.
G1-201LBX1	Bariatric steel folding commode	650 lbs. (294.8 kg)	16"-22"	16.5"	19"	21"	1 ea
MDS89668XW	Bariatric steel fixed drop arm bench seat commode	850 lbs. (385.5 kg)	21"	21"	26.5"	29"	1 ea
MDS89668XW1K	Bariatric steel fixed drop arm bench seat commode	1,000 lbs. (453.6 kg)	19.75"-23.75"	21"	26.5"	29"	1 ea
G1-400BX2	Bariatric steel fixed bench seat commode	650 lbs. (294.8 kg)	17 "-20.5"	18.5"	24"	29.5"	2 ea



G1-202ELX4



G1-201LBX1



MDS89668XW1K



G1-400BX2

## Commode bucket

- Replacement commode bucket includes lid and metal handle

Item No.	Description	Pkg.
MDS80306B	Grey	6 cs
MDS80306EB	Grey bucket, no lid, no handle	18 cs



MDS80306B

### Shower chairs with back

- Rustproof aluminum shower chairs for patients who have difficulty standing in the shower; backrest included for added stability, security and comfort

Item No.	Description	Wt. capacity	Seat height	Seat depth	Pkg.
MDS89745RAH	White shower chair with arms	350 lbs. (158.8 kg)	16"-21"	13.5	1 ea
MDS89745A	White assembled shower chair	250 lbs. (113.4 kg)	13.75"-20.75"	12.25	2 ea
MDS89745ELMBH	Tan shower chair with arms, Microban*	350 lbs. (158.8 kg)	15"-21"	15	1 ea
MDS89745KDMBH	Light blue shower chair, Microban*	300 lbs. (136.1 kg)	13.75"-20.75"	12.25	1 ea
G2-100KHX1	Grey shower chair with arms	350 lbs. (158.8 kg)	17.25"-22.5"	13.5"	1 ea
G2-101KX1	Grey tool-free shower chair	400 lbs. (181.4 kg)	13.75"-20.5"	12	1 ea
G2-101KRX1	Grey retail shower chair, tool-free	400 lbs. (181.4 kg)	13.75"-20.5"	12	1 ea
G2-102KX1	Grey knockdown shower chair	400 lbs. (181.4 kg)	15.25"-19"	12	1 ea
G2-102KX4	Grey knockdown shower chair	400 lbs. (181.4 kg)	15.25"-19"	12	4 ea

\*Antimicrobial properties are built in to protect the product. The product does not protect users or others against bacteria, viruses, germs or other disease-causing organisms



MDS89745RAH



G2-102KX1

### Shower chairs without back

- Rust-proof aluminum shower chairs for patients who have difficulty standing in the shower; no backrest for easy maneuverability in tight spaces

Item No.	Description	Wt. capacity	Seat height	Seat depth	Pkg.
MDS89740RWAH	White shower chair with arms	350 lbs. (158.8 kg)	16"-21"	13.5"	1 ea
MDS89740A	White assembled shower chair	250 lbs. (113.4 kg)	13.75"-20.75"	12"	2 ea
G2-201KRX1	Grey retail shower chair, tool-free	400 lbs. (181.4 kg)	13.75"-20.5"	12"	1 ea
G2-200KHX1	Grey shower chair with arms	350 lbs. (158.8 kg)	17.25"-22.5"	13.5"	1 ea
G2-202KX1	Grey shower chair, knockdown	400 lbs. (181.4 kg)	15"-19"	12"	1 ea
G2-300KRX1	Grey bath stool	300 lbs. (136.1 kg)	13.75"-20.5"	12.5"	1 ea
G2-300KRPX1	Grey padded bath stool with swivel	300 lbs. (136.1 kg)	16"-22"	13.75"	1 ea



MDS89740RWAH



G2-201KRX1



G2-300KRPX1



### Bariatric shower chairs

- Heavy duty reinforced frames for high weight capacity needs

Item No.	Description	Wt. capacity	Seat height	Seat depth	Pkg.
MDS89745AXW	Bench with back	550 lbs. (249.5 kg)	16.5"-20.5"	20.5"	1 ea
MDS89740AXW	Bench without back	550 lbs. (249.5 kg)	14"-17"	12"	1 ea
G2-102BX1	White bench with back	550 lbs. (249.5 kg)	15"-20"	12"	1 ea
G2-202BX1	White bench without back	550 lbs. (249.5 kg)	15"-20"	12"	1 ea

### Transfer benches

- Allows patients to enter the bath tub from a seated position; padded and bariatric options available to meet various needs

Item No.	Description	Wt. capacity	Seat height	Seat depth	Pkg.
G3-100KBX1	Grey transfer bench, Brown box	400 lbs. (181.4 kg)	18"-22.5"	16.25"	1 ea
G3-100KBX2	Grey transfer bench, Brown box	400 lbs. (181.4 kg)	18"-22.5"	16.25"	2 cs
G3-100KBRX1	Grey transfer bench, retail box	400 lbs. (181.4 kg)	18"-22.5"	16.25"	1 ea
G3-200KBX1	Grey deluxe transfer bench, Brown box	500 lbs. (226.8 kg)	19.25"-23.5"	16.5"	1 ea
MDS86952XW	White premium bariatric	550 lbs. (249.5 kg)	16"-20"	16"	1 ea



### Elevated toilet seats

- Adds up to 5" in height to toilet seats for patients with difficulty sitting and standing; available with or without arms

Item No.	Description	Wt. capacity	Pkg.
MDS80314	White elevated locking toilet seat, without arms	350 lbs. (158.8 kg)	3 cs
MDS80314H	White elevated locking toilet seat, without arms	350 lbs. (158.8 kg)	1 ea
MDS80315EL	White elongated elevated locking toilet seat with arms	350 lbs. (158.8 kg)	3 cs
MDS80316	White elevated locking toilet seat, with arms	350 lbs. (158.8 kg)	3 cs
MDS80316H	White elevated locking toilet seat, with arms	350 lbs. (158.8 kg)	1 ea
MDS80316MB	Blue elevated locking toilet seat, with arms, with Microban*	350 lbs. (158.8 kg)	3 cs
G30250	White toilet seat riser, without arms, without locking, 5"	250 lbs. (113.4 kg)	3 cs
G30250H	White toilet seat riser, without arms, without locking, 5"	250 lbs. (113.4 kg)	1 ea
G30260	White toilet seat riser, without arms, with locking, 5"	350 lbs. (158.8 kg)	3 cs
G30270A	White toilet seat riser, with arms, with locking, 5"	350 lbs. (158.8 kg)	3 cs
G30270AH	White toilet seat riser, with arms, with locking, 5"	350 lbs. (158.8 kg)	1 ea

\*Antimicrobial properties are built in to protect the product. The product does not protect users or others against bacteria, viruses, germs or other disease-causing organisms



G30250



G30270A



G30260



MDS80315EL



MDS80316MB

## White grab bars

- Enamel coated steel grab bars mount to the wall

Item No.	Description	Dimensions	Wt. capacity	Pkg.
MDS86012	White enamel grab bar	12" wide	250 lbs. (113.4 kg)	3 cs
MDS86012H	White enamel grab bar	12" wide	250 lbs. (113.4 kg)	1 ea
G87520M	White enamel grab bar	12" wide	300 lbs. (136.1 kg)	3 cs
MDS86016	White enamel grab bar	16" wide	250 lbs. (113.4 kg)	3 cs
MDS86016H	White enamel grab bar	16" wide	250 lbs. (113.4 kg)	1 ea
G87530M	White enamel grab bar	16" wide	300 lbs. (136.1 kg)	3 cs
MDS86018	White enamel grab bar	18" wide	250 lbs. (113.4 kg)	3 cs
MDS86018H	White enamel grab bar	18" wide	250 lbs. (113.4 kg)	1 ea
G87540M	White enamel grab bar	18" wide	300 lbs. (136.1 kg)	3 cs
MDS86024	White enamel grab bar	24" wide	250 lbs. (113.4 kg)	3 cs
G87550M	White enamel grab bar	24" wide	300 lbs. (136.1 kg)	3 cs
MDS86032	White enamel grab bar	32" wide	250 lbs. (113.4 kg)	3 cs
G87560M	White enamel grab bar	32" wide	300 lbs. (136.1 kg)	3 cs
G87570M	White enamel grab bar	36" wide	300 lbs. (136.1 kg)	3 cs



## Chrome grab bars

- Chrome plated steel grab bars mount to the wall

Item No.	Description	Dimensions	Wt. capacity	Pkg.
MDS86012CHR	Knurled chrome grab bar	12" wide	250 lbs. (113.4 kg)	3 cs
MDS86012CHRH	Knurled chrome grab bar	12" wide	250 lbs. (113.4 kg)	1 ea
G81020KEM	Knurled chrome grab bar	12" wide	300 lbs. (136.1 kg)	3 cs
MDS86016CHR	Knurled chrome grab bar	16" wide	250 lbs. (113.4 kg)	3 cs
MDS86016CHRH	Knurled chrome grab bar	16" wide	250 lbs. (113.4 kg)	1 ea
G81030KEM	Knurled chrome grab bar	16" wide	300 lbs. (136.1 kg)	3 cs
MDS86018CHR	Knurled chrome grab bar	18" wide	250 lbs. (113.4 kg)	3 cs
MDS86018CHRH	Knurled chrome grab bar	18" wide	250 lbs. (113.4 kg)	1 ea
G81040KEM	Knurled chrome grab bar	18" wide	300 lbs. (136.1 kg)	3 cs
MDS86024CHR	Knurled chrome grab bar	24" wide	250 lbs. (113.4 kg)	3 cs
MDS86024CHRH	Knurled chrome grab bar	24" wide	250 lbs. (113.4 kg)	1 ea
MDS86032CHR	Knurled chrome grab bar	32" wide	250 lbs. (113.4 kg)	3 cs
MDS86032CHRH	Knurled chrome grab bar	32" wide	250 lbs. (113.4 kg)	1 ea
G81050KEM	Knurled chrome grab bar	24" wide	300 lbs. (136.1 kg)	3 cs
G81060KEM	Knurled chrome grab bar	32" wide	300 lbs. (136.1 kg)	3 cs
G81070KEM	Knurled chrome grab bar	36" wide	300 lbs. (136.1 kg)	3 cs





### Suction grab bars

- Attach to the wall without screws

Item No.	Description	Wt. capacity	Pkg.
MDS86012SCLI	White with locking indicator, 12"L	250 lbs. (113.4 kg)	3 cs
MDS86016SCLI	White with locking indicator, 16"L	250 lbs. (113.4 kg)	3 cs
MDS86019SCLI	White, pivot grab bar, 19"L	250 lbs. (113.4 kg)	3 cs



MDS86012SCLI

### Tub grab bars

- Securely mount to the tub wall

Item No.	Description	Wt. capacity	Pkg.
MDS86321W	White, locks to side of tub	250 lbs. (113.4 kg)	2 cs
MDS86321WH	White, locks to side of tub	250 lbs. (113.4 kg)	1 ea
G98006	White, bi-level grab bar	250 lbs. (113.4 kg)	3 cs



MDS86321W

### Cast protectors

- Help keep casts dry when bathing; available in adult and child sizes

Item No.	Description	Pkg.
CUR200ALL	Clear, adult leg	6 cs
CUR200ALLH	Clear, adult leg	1 ea
CUR200AAA	Clear, adult arm	6 cs
CUR200AAAH	Clear, adult arm	1 ea
CUR100KLL	Clear, child leg	6 cs
CUR100KLLH	Clear, child leg	1 ea
CUR100KAA	Clear, child arm	6 cs
CUR100KAAH	Clear, child arm	1 ea



CUR200ALL



CUR200AAA

## Shower heads

- Easy to install white handheld shower heads with diverter valve

Item No.	Description	Pkg.
MDSHHSOWER	6.5"-ft. long hose	6 cs
MDSHHSOWERH	6.5"-ft. long hose	1 ea
G98035	Handheld shower hose, 7.2"-ft. long hose	4 cs
G98035H	Handheld shower hose, 7.2"-ft. long hose	1 ea

## Bath mats

- Suction cup rubber bath mats with Microban® help prevent slipping

Item No.	Description	Pkg.
MDSBATHMAT	White bath mat with Microban*	4 cs
MDSBATHMATH	White bath mat with Microban*	1 ea

\*Antimicrobial properties are built in to protect the product. The product does not protect users or others against bacteria, viruses, germs or other disease-causing organisms



MDSHHSOWER



MDSBATHMATH

# Wheelchairs

## K1 Guardian wheelchairs

- Durable hammertone frame; available in nylon or vinyl upholstery
- Hemi-height adjustable
- Accommodates 2" seat extender

### Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length swing back	Swing-away	K1166N22S	Black	300 lbs. (136.1 kg)	16" x 16"
		K1186N22S	Black	300 lbs. (136.1 kg)	18" x 16"
		K1206N22S	Black	300 lbs. (136.1 kg)	20" x 16"
	Elevating	K1166N22E	Black	300 lbs. (136.1 kg)	16" x 16"
		K1186N22E	Black	300 lbs. (136.1 kg)	18" x 16"
		K1206N22E	Black	300 lbs. (136.1 kg)	20" x 16"
Full-length removable	Swing-away	K1186N13S	Black	300 lbs. (136.1 kg)	18" x 16"
	Elevating	K1186N13E	Black	300 lbs. (136.1 kg)	18" x 16"

### Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length swing back	Swing-away	K1166V22S	Black	300 lbs. (136.1 kg)	16" x 16"
		K1186V22S	Black	300 lbs. (136.1 kg)	18" x 16"
		K1206V22S	Black	300 lbs. (136.1 kg)	20" x 16"
	Elevating	K1166V22E	Black	300 lbs. (136.1 kg)	16" x 16"
		K1186V22E	Black	300 lbs. (136.1 kg)	18" x 16"
		K1206V22E	Black	300 lbs. (136.1 kg)	20" x 16"
Full-length removable	Swing-away	K1186V11S	Black	300 lbs. (136.1 kg)	18" x 16"
	Elevating	K1186V11E	Black	300 lbs. (136.1 kg)	18" x 16"

Swing-away arm rest

Accommodates 2" seat depth extenders



K1166V22S

## K1 Basic wheelchairs

- Durable hammertone frame
- Hemi height adjustable

### Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length removable	Swing-away	MDS806250NEE	Black	300 lbs. (136.1 kg)	16" x 16"
		MDS806250EE	Black	300 lbs. (136.1 kg)	18" x 16"
		MDS806250EET	Black	300 lbs. (136.1 kg)	18" x 16"
	Elevating	MDS806400EE	Black	300 lbs. (136.1 kg)	20" x 16"
		MDS806300NEE	Black	300 lbs. (136.1 kg)	16" x 16"
		MDS806300EE	Black	300 lbs. (136.1 kg)	18" x 16"
Full-length permanent	Elevating	MDS806450EE	Black	300 lbs. (136.1 kg)	20" x 16"
	Swing-away	MDS806150EE	Black	300 lbs. (136.1 kg)	18" x 16"
	Elevating	MDS806200EE	Black	300 lbs. (136.1 kg)	18" x 16"



MDS806300NEE

## K2 Guardian wheelchairs

- Durable hammertone frame
- Ships in hemi height; includes seat belt, anti tippers and brake extenders

### Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Full-length removable and height adjustable	Elevating	K2166N42E1	Black	300 lbs. (136.1 kg)	16" x 16"
		K2186N42E1	Black	300 lbs. (136.1 kg)	18" x 16"
		K2206N42E1	Black	300 lbs. (136.1 kg)	20" x 16"

Accommodates  
2" seat  
extenders



K2166N42E1

## K2 Basic wheelchairs

- Durable tig-welded frame with silver vein finish
- Hemi height adjustable axles raise and lower seat 2"
- Easy-to-clean vinyl upholstery

### Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, removable	Swing-away	MDS806250NEV	Black	300 lbs. (136.1 kg)	16" x 16"
		MDS806250EV	Black	300 lbs. (136.1 kg)	18" x 16"
		MDS806400EV	Black	300 lbs. (136.1 kg)	20" x 16"
	Elevating	MDS806300NEV	Black	300 lbs. (136.1 kg)	16" x 16"
		MDS806300EV	Black	300 lbs. (136.1 kg)	18" x 16"
		MDS806450EV	Black	300 lbs. (136.1 kg)	20" x 16"
Full-length permanent	Swing-away	MDS806150EV	Black	300 lbs. (136.1 kg)	18" x 16"
	Elevating	MDS806200EV	Black	300 lbs. (136.1 kg)	18" x 16"



MDS806300EV

### K3 Basic wheelchairs

- Durable hammertone frame
- Hemi height adjustable axles

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, removable	Swing-away	MDS806600NE	Black	300 lbs. (136.1 kg)	16" x 16"
		MDS806600E	Black	300 lbs. (136.1 kg)	18" x 16"
		MDS806660E	Black	300 lbs. (136.1 kg)	20" x 16"
	Elevating	MDS806650NE	Black	300 lbs. (136.1 kg)	16" x 16"
		MDS806650E	Black	300 lbs. (136.1 kg)	18" x 16"
		MDS806665E	Black	300 lbs. (136.1 kg)	20" x 16"



MDS806600E

### K3 Plus wheelchairs

- Durable hammertone frame
- Hemi height adjustable axles

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, swing back and height adjustable	Swing-away	MDS806600NEPL	Black	300 lbs. (136.1 kg)	16" x 16"
		MDS806600EPL	Black	300 lbs. (136.1 kg)	18" x 16"
		MDS806660EPL	Black	300 lbs. (136.1 kg)	20" x 16"
	Elevating	MDS806650NEPL	Black	300 lbs. (136.1 kg)	16" x 16"
		MDS806650EPL	Black	300 lbs. (136.1 kg)	18" x 16"
		MDS806665EPL	Black	300 lbs. (136.1 kg)	20" x 16"



MDS806650EPL

### K4 lightweight wheelchairs

- Durable hammertone frame

#### Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, swing back	Swing-away	MDS806500N	Black	300 lbs. (136.1 kg)	16" x 16"
		MDS806500	Black	300 lbs. (136.1 kg)	18" x 16"
		MDS806560	Black	350 lbs. (158.8 kg)	20" x 16"
	Elevating	MDS806570	Black	350 lbs. (158.8 kg)	22" x 16"
		MDS806550N	Black	300 lbs. (136.1 kg)	16" x 16"
		MDS806550	Black	300 lbs. (136.1 kg)	18" x 16"
		MDS806565	Black	350 lbs. (158.8 kg)	20" x 16"
		MDS806575	Black	350 lbs. (158.8 kg)	22" x 16"
		MDS806500FLA	Black	300 lbs. (136.1 kg)	18" x 16"
Full-length removable	Swing-away	MDS806500FLA	Black	300 lbs. (136.1 kg)	18" x 16"
	Elevating	MDS806550FLA	Black	300 lbs. (136.1 kg)	18" x 16"



MDS806560

## K4 Basic wheelchairs

- Durable grey frame
- Hemi height adjustable

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, swing back	Swing-away	MDS806500NE	Silver	300 lbs. (136.1 kg)	16" x 16"
		MDS806500E	Silver	300 lbs. (136.1 kg)	18" x 16"
		MDS806500NE18	Silver	300 lbs. (136.1 kg)	16" x 18"
	Elevating	MDS806560E	Silver	300 lbs. (136.1 kg)	20" x 18"
		MDS806550NE	Silver	300 lbs. (136.1 kg)	16" x 16"
		MDS806550E	Silver	300 lbs. (136.1 kg)	18" x 16"
		MDS806500E18	Silver	300 lbs. (136.1 kg)	18" x 18"
		MDS806565E	Silver	300 lbs. (136.1 kg)	20" x 18"



MDS806500E

## K4 Plus wheelchairs

- Durable hammertone frame
- Includes antitippers

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, swing-back and height adjustable	Swing-away	MDS806500PLUS	Black	300 lbs. (136.1 kg)	18" x 16"
	Elevating	MDS806560PLUS	Black	300 lbs. (136.1 kg)	20" x 18"
		MDS806550PLUS	Black	300 lbs. (136.1 kg)	18" x 16"



Includes anti tippers



MDS806550PLUS

## Guardian recliners

- Infinite position reclining up to 180°
- Durable hammertone frame

## Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, removable	Elevating	KR168N21E	Black	300 lbs. (136.1 kg)	16" x 18"
		KR188N21E	Black	300 lbs. (136.1 kg)	18" x 18"
		KR208N21E	Black	300 lbs. (136.1 kg)	20" x 18"
Full-length removable	Elevating	KR188N11E	Black	300 lbs. (136.1 kg)	18" x 18"



KR168N21E

### Reclining wheelchairs

- Infinite position reclining from 90°–140°
- Durable carbon steel frame
- Hemi-height adjustable

#### Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, removable	Elevating	MDS808350	Chrome	300 lbs. (136.1 kg)	16" x 17"
		MDS808450	Chrome	300 lbs. (136.1 kg)	18" x 17"
		MDS808550	Chrome	300 lbs. (136.1 kg)	20" x 17"
		MDS808650	Chrome	300 lbs. (136.1 kg)	22" x 17"



MDS808550

### Hybrid wheelchairs

- Combination wheelchair and transport chair
- Quick release large wheel

#### Nylon

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, removable	Swing-away	MDS806250NH2	Blue	300 lbs. (136.1 kg)	16" x 16"
	Elevating	MDS806250H2	Blue	300 lbs. (136.1 kg)	18" x 16"
		MDS806300H2	Blue	300 lbs. (136.1 kg)	18" x 16"



MDS806250H2

### Excel wheelchairs

- Durable carbon steel frame
- Hemi-height adjustable

#### Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, removable	Swing-away	MDS806250N	Chrome	300 lbs. (136.1 kg)	16" x 16"
		MDS806250D	Chrome	300 lbs. (136.1 kg)	18" x 16"
		MDS806250NVY	Chrome	300 lbs. (136.1 kg)	18" x 16"
		MDS806400	Chrome	300 lbs. (136.1 kg)	20" x 16"
	Elevating	MDS806300N	Chrome	300 lbs. (136.1 kg)	16" x 16"
		MDS806300D	Chrome	300 lbs. (136.1 kg)	18" x 16"
		MDS806300NVY	Chrome	300 lbs. (136.1 kg)	18" x 16"
		MDS806450	Chrome	300 lbs. (136.1 kg)	20" x 16"
Full-length permanent	Swing-away	MDS806150D	Chrome	300 lbs. (136.1 kg)	18" x 16"
	Elevating	MDS806200N	Chrome	300 lbs. (136.1 kg)	16" x 16"
Full-length removable	Swing-away	MDS806100D	Chrome	300 lbs. (136.1 kg)	18" x 16"
		MDS806400FLA	Chrome	300 lbs. (136.1 kg)	20" x 16"
	Elevating	MDS806300DFLA	Chrome	300 lbs. (136.1 kg)	18" x 16"



MDS806300D

### Excel extra-wide wheelchairs

- Durable carbon steel frame
- Hemi-height adjustable

#### Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, removable	Swing-away	MDS806700	Chrome	500 lbs. (226.8 kg)	20" x 18"
		MDS806800	Chrome	500 lbs. (226.8 kg)	22" x 18"
		MDS806900	Chrome	500 lbs. (226.8 kg)	24" x 18"
	Elevating	MDS806750	Chrome	500 lbs. (226.8 kg)	20" x 18"
		MDS806850	Chrome	500 lbs. (226.8 kg)	22" x 18"
		MDS806950	Chrome	500 lbs. (226.8 kg)	24" x 18"
Full-length removable	Swing-away	MDS806700FLA	Chrome	500 lbs. (226.8 kg)	20" x 18"
		MDS806900FLA	Chrome	500 lbs. (226.8 kg)	24" x 18"
		MDS806750FLA	Chrome	500 lbs. (226.8 kg)	20" x 18"
	Elevating	MDS806850FLA	Chrome	500 lbs. (226.8 kg)	22" x 18"
		MDS806950FLA	Chrome	500 lbs. (226.8 kg)	24" x 18"



MDS806900

### Shuttle bariatric wheelchairs

- Durable carbon steel frame
- Hemi-height adjustable

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, removable	Elevating	MDS809650	Black	850 lbs. (385.5 kg)	24" x 20"
		MDS809750	Black	850 lbs. (385.5 kg)	26" x 20"
		MDS809850	Black	850 lbs. (385.5 kg)	28" x 20"
		MDS809950	Black	850 lbs. (385.5 kg)	30" x 20"



MDS809750

### Pediatric wheelchair

- Adjustable telescoping handles.
- Hemi-height adjustable

#### Vinyl

Arm style	Leg style	Item No.	Frame color	Wt. capacity	Seat dimensions
Desk-length, swing-back	Elevating	MDS806140PD	Black	250 lbs. (113.4 kg)	14" x 12"



MDS806140PD



## Wheelchair accessories

Item No.	Description	Pkg.
MDS85190	Holder, oxygen tank, IV pole, wheelchair	1 ea
WCA6100	Holder, oxygen tank, for wheelchair ADL	1 ea
MDS85189	Device, anti-tip, wheelchair, rear	1 pr
STMSM22	Anti-rollback device	1 ea
MDS85189H	Bag, leg rest or carry, all for wheelchairs	1 pr
MDS85181U	Universal oxygen tank holder for wheelchairs	1 ea
MDS85189H	Anti-tip device, K1, rear, hemi	1 ea
MDS85183	Wheelchair IV pole, telescoping	1 ea
MDS85189WK4	Anti-tip rear wheel for K4 and recliner	1 pr
WCABELT	Safety belt for wheelchair	1 ea
MDS85157PR	Black handle extension wheel lock 5/8 x 8	1 pr
MDS85197A	Adjustable anti-theft device for wheelchair	1 ea
MDSK41RMDRIVE	One arm drive attachment for Excel K	1 ea
WCA806991	Nut, tinnerman for wheelchair, leg rigging	1 pr
WCA6200	Cup holder, ADL	1 ea
MDS85196	Anti-fold device	1 ea
MDS85195	Kit, super hemi, with standard wheel	1 ea
MDS851965SH	Anti-fold device for shuttle wheelchair	1 ea
MDS85197SH	Anti-theft, overhead for shuttle wheelchair	1 ea
STMSM22A	Anti-rollback device with alarm	1 ea
MDS85180	5-in-1 IV/O2/anti-theft device	1 ea
STMSM22W	Anti-rollback device	1 ea
MDS85189SH	Anti-tip, bariatric, for shuttle wheelchair	1 pr
MDS85189K	Anti-tip Kidz chair	1 pr
WCA806922EXT	K4 seat extension. 18" x 18" 3000 series	1 ea
MDS806600EXT	2 exit kit for K3 Basic Plus, 18"	1 ea
WCAA0060	K1/K2 nylon seat extender, 18"	1 ea
WCAA0061	K1/K2 nylon seat extender, 20"	1 ea



# Transport chairs

## Basic steel transport chairs with 8" wheels

- Durable steel frame with black upholstery
- Includes detachable footrests and a seat belt for safety
- Back folds down for convenient transporting

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	Pkg.
MDS808200B	Black	250 lbs. (113.4 kg)	Swing-away	Full-length permanent	19"	1 ea
MDS808200	Chrome	300 lbs. (136.1 kg)	Swing-away	Full-length permanent	19"	1 ea



MDS808200B

## Aluminum transport chairs with 8" wheels

- Powder coated aluminum frame with black upholstery
- Includes seat belt for safety and padded arms

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	Pkg.
MDS808200ABE	Blue	300 lbs. (136.1 kg)	Swing-away	Full-length permanent	19"	1 ea
MDS808200ARE	Red	300 lbs. (136.1 kg)	Swing-away	Full-length permanent	19"	1 ea



MDS808200ABE

## Aluminum transport chairs with 12" wheels

- Powder coated aluminum frame with black upholstery
- Includes seat belt for safety and padded arms
- Larger back tire for maneuvering through obstacles
- Hand brakes on handles

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	Pkg.
MDS808210ABE	Blue	300 lbs. (136.1 kg)	Swing-away	Full-length permanent	19"	1 ea
MDS808210ARE	Red	300 lbs. (136.1 kg)	Swing-away	Full-length permanent	19"	1 ea



MDS808210ARE

### Excel® Translator combination rollator/transport chairs

- Durable steel frame with black upholstery
- Includes easy-to-adjust, push-button, swing-away footrests
- Footrests lock to sides of rollator when not in use
- Height-adjustable push grips and comfortable hand brake



MDS808200TR

Item No.	Description	Pkg.
MDS808200TRR	Red	1 ea
MDS808200TR	Blue	1 ea

### Ultralight transport chairs

- Lightweight aluminum frame with black upholstery
- Carrying handles and compact design allow for easy and convenient storage and transport



MDS808200F3B

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	Pkg.
MDS808200F3R	Red	300 lbs. (136.1 kg)	Swing-away	Desk-length permanent	19"	1 ea
MDS808200F3B	Blue	300 lbs. (136.1 kg)	Swing-away	Desk-length permanent	19"	1 ea
MDS808200F3BK	Black	300 lbs. (136.1 kg)	Swing-away	Desk-length permanent	19"	1 ea
MDS808200F3S	Silver	300 lbs. (136.1 kg)	Swing-away	Desk-length permanent	19"	1 ea

### Bariatric transport chairs

- Durable frame with black upholstery, seat belt, handbrakes and wheel lock for better control



MDS80822B

Item No.	Description	Wt. capacity	Leg style	Arm style	Seat width	Pkg.
MDS80822B	Bariatric steel, Red	500 lbs. (226.8 kg)	Swing-away	Desk-length	22"	1 ea
MDS808200BAR	Bariatric, lightweight, aluminum, Red	400 lbs. (181.4 kg)	Elevating	Elevating	22"	1 ea

# Rollators

## Basic steel rollators

- Durable powder-coated steel frame
- Storage bag under the seat is included

Item No.	Description	Wt. capacity	Seat height	Wheel size	Height range	Hand grip width	Pkg.
MDS86850ES	Burgundy	350 lbs. (158.8 kg)	20.5"	6"	31"-35"	17.5"	1 ea
MDS86850EBS	Blue	350 lbs. (158.8 kg)	20.5"	6"	31"-35"	17.5"	1 ea
MDS86850ESKD	Burgundy	350 lbs. (158.8 kg)	20.5"	6"	31"-35"	17"	1 ea
MDS86830EBLS	Black	300 lbs. (136.1 kg)	20.5"	6"	31"-35"	17.5"	1 ea
MDS86850ES8	Burgundy	350 lbs. (158.8 kg)	20.5"	8"	31"-35"	17.5"	1 ea
MDS86850EBS8	Blue	350 lbs. (158.8 kg)	20.5"	8"	31"-35"	17.5"	1 ea
MDS86830EBLS8	Black	300 lbs. (136.1 kg)	20.5"	8"	31"-35"	17.5"	1 ea
MDS86840EBS	Blue	300 lbs. (136.1 kg)	20.5"	6"	31"-36.25"	19.75"	1 ea
MDS86840EBS8	Blue	300 lbs. (136.1 kg)	21"	8"	31.25"-36.25"	19.75"	1 ea



MDS86850ES8

## Basic aluminum rollators

- Durable aluminum frame
- Storage bag under the seat is included

Item No.	Description	Wt. capacity	Seat height	Wheel size	Height range	Hand grip width	Pkg.
MDS86850EB	Blue	250 lbs. (113.4 kg)	20.5"	6"	31"-35"	17.5"	1 ea
MDS86850E	Burgundy	250 lbs. (113.4 kg)	20.5"	6"	31"-35"	17.5"	1 ea
MDS86850EP	Purple	250 lbs. (113.4 kg)	20.5"	6"	31"-35"	17.5"	1 ea
MDS86850EG	Green	250 lbs. (113.4 kg)	20.5"	6"	31"-35"	17.5"	1 ea
MDS86830EBL	Black	300 lbs. (136.1 kg)	20.5"	6"	31"-35"	17.5"	1 ea
MDS86850EJR	Burgundy	250 lbs. (113.4 kg)	20.5"	6"	29"-33"	17.5"	1 ea
MDS86830EBL8	Black	350 lbs. (158.8 kg)	20.5"	8"	31"-35"	17.5"	1 ea



MDS86850EB

## Superlight rollators with 6" wheels

- Lightweight aluminum frame
- Weighs only 11 lbs. making it easy to transport or store when not in use

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	Pkg.
MDS86825SLR	Burgundy	250 lbs. (113.4 kg)	20.5"	29.5"-36"	17"	1 ea
MDS86825SLBLK	Blue	250 lbs. (113.4 kg)	20.5"	29.5"-36"	17"	1 ea
MDS86825SLB	Purple	250 lbs. (113.4 kg)	20.5"	29.5"-36"	17"	1 ea
MDS86825SLRS*	Red	350 lbs. (158.8 kg)	18"-23"	29.5"-36"	17"	1 ea



MDS86825SLBLK

Aluminum options weigh only 11 lbs!

\*MDS86825SLRS made of steel

### Momentum rollators with 6" wheels

- Lightweight steel frame with large padded seat
- Seat flips up to enable user to walk inside
- All-terrain wheels make it easier to get around
- Product weighs 18.8 lbs.

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	Pkg.
MDS86870R	Red	250 lbs. (113.4 kg)	20.5"	34.25"–41.5"	18"	1 ea
MDS86870W	White	250 lbs. (113.4 kg)	20.5"	34.25"–41.5"	18"	1 ea
MDS86870G	Grey	250 lbs. (113.4 kg)	20.5"	34.25"–41.5"	18"	1 ea

### Empower rollators with 8" wheels

- Aluminum frame with memory foam seat
- Padded backrest flips forward and backward enabling the user to sit comfortably on either side
- Seat flips up to enable user to walk inside

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	Pkg.
MDS86845BM	Blue	300 lbs. (136.1 kg)	20.5"	34.25"–41.5"	18"	1 ea
MDS86845BLKM	Black	300 lbs. (136.1 kg)	20.5"	34.25"–41.5"	18"	1 ea
MDS86845SM	Silver	300 lbs. (136.1 kg)	20.5"	34.25"–41.5"	18"	1 ea

### Army rollator with 6" wheels

- Durable powder-coated steel frame
- Black cooler bag with Army Star logo

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	Pkg.
MDS86850ESKDA	Black	350 lbs. (158.8 kg)	20.5"	31"–35"	17.5"	1 ea



Strong in mind, body, and soul. The U.S. Army supports our nation and our values every day. Now, you can support Army service members past and present with U.S. Army licensed gear from Medline. Learn more about what it means to be a U.S. Army soldier at [www.goarmy.com](http://www.goarmy.com).



### Bariatric rollators with 8" wheels

- Heavy-duty bariatric rollators

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	Pkg.
MDS86800XWS	Steel, Burgundy	500 lbs. (226.8 kg)	20.5"	31.5"-37.5"	23.5"	1 ea
MDS86800XW	Aluminum, Burgundy	400 lbs. (181.4 kg)	20.5"	31.5"-37.5"	23.5"	1 ea

### Guardian rollators with 8" wheels

- Aluminum frame
- Padded seat and curved padded backrest provide enhanced support and comfort

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	Pkg.
G07887B	Blue	300 lbs. (136.1 kg)	20"	31"-36.5"	17.5"	1 ea
G07887R	Red	300 lbs. (136.1 kg)	20"	31"-36.5"	17.5"	1 ea

### Steel rollators with basket

- Easily removable wired basket and a flip up backrest for easy storage.

Item No.	Description	Wt. capacity	Seat height	Height range	Hand grip width	Pkg.
MDS86860ERS8	Red	300 lbs. (136.1 kg)	21.25"	30.5"-36"	18.5"	1 ea
MDS86860EBS8	Blue	300 lbs. (136.1 kg)	21.25"	30.5"-36"	18.5"	1 ea



Flip-up  
back rest

## byACRE carbon fiber rollators

- World's lightest rollators at 10.6 lbs.
- Carbon fiber frame is 5X stronger than steel
- Ergonomic handles have integrated brake lines and are height adjustable in 1" increments
- All byACRE rollators include organizer bag

### Compact

Item No.	Description	Wt. capacity	Seat height	Height range	Seat width	Pkg.
BYA100SMBK	Black	285 lbs. (129.3 kg)	22"	31"-36"	16"	1 ea
BYA100SMR	Red	285 lbs. (129.3 kg)	22"	31"-36"	16"	1 ea
BYA100SMW	White	285 lbs. (129.3 kg)	22"	31"-36"	16"	1 ea

### Regular

Item No.	Description	Wt. capacity	Seat height	Height range	Seat width	Pkg.
BYA100MDBK	Black	285 lbs. (129.3 kg)	24"	31"-36"	16.5"	1 ea
BYA100MDR	Red	285 lbs. (129.3 kg)	24"	31"-36"	16.5"	1 ea
BYA100MDW	White	285 lbs. (129.3 kg)	24"	31"-36"	16.5"	1 ea

### Comfort

Item No.	Description	Wt. capacity	Seat height	Height range	Seat width	Pkg.
BYA100LGBK	Black	285 lbs. (129.3 kg)	24"	31"-36"	18.5"	1 ea
BYA100LGR	Red	285 lbs. (129.3 kg)	24"	31"-36"	18.5"	1 ea
BYA100LGW	White	285 lbs. (129.3 kg)	24"	31"-36"	18.5"	1 ea



BYA100SMBK



BYA100MDW



BYA100LGR

Folds for easy portability



**byACRE carbon fiber rollator accessories**

Item No.	Description	Pkg.
BYAGRBG	Grocery bag, Grey	1 ea
BYAWKNBG	Weekend bag, Grey	1 ea
BYATVLBG	Travel bag, Grey	1 ea
BYACNHR	Cane holder, right	1 ea
BYACNHL	Cane holder, left	1 ea
BYAOBGR	Organizer bag, small/medium, Grey	1 ea
BYAOBGL	Organizer bag, large, Grey	1 ea
BYASMBR	Backrest, small/medium	1 ea
BYALBR	Backrest, large	1 ea



BYAGRBG



BYAWKNBG



BYATVLBG



BYAOBGR



BYASMBR





# Walkers

## Two-button folding walkers\*

- Durable walker with two-button folding system allows easy movement through spaces

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	Pkg.
MDS864104	Aluminum, no wheels	300 lbs. (136.1 kg)	PVC	32"-39"	18"	4 cs
MDS86410W4	Aluminum, 3" wheels	300 lbs. (136.1 kg)	PVC	32"-38"	17.5"	4 cs
MDS86410W54	Aluminum, 5" wheels	300 lbs. (136.1 kg)	PVC	32"-38"	17"	4 cs
MDS86410W54B	Steel, 5" wheels	400 lbs. (181.4 kg)	PVC	32.5"-39.5"	17"	4 cs

## Youth two-button folding walkers

- Youth size durable walker with two-button folding system allows easy movement through spaces

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	Pkg.
MDS86410J4	Aluminum, no wheels	300 lbs. (136.1 kg)	PVC	25"-32"	17.5"	4 cs
MDS86410JW54	Aluminum, 5" wheels	300 lbs. (136.1 kg)	PVC	25"-32"	17.5"	4 cs
MDS86410JW54B	Steel, 5" wheels	400 lbs. (181.4 kg)	PVC	25"-32"	17.5"	4 cs



MDS864104



MDS86410W54



MDS86410JW54

\*All MDS walkers come unassembled

## Bariatric folding walkers

- Extra-wide frames and higher weight capacities provide a comfortable fit for larger patients
- Includes two-button folding system
- All bariatric folding walkers come assembled

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	Pkg.
MDS86410XW	Aluminum, no wheels	500 lbs. (226.8 kg)	PVC	31"-41"	20.5"	1 ea
MDS86410XWB	Steel, no wheels	650 lbs. (294.8 kg)	PVC	31"-41"	20"	1 ea
MDS86410XWW	Aluminum, 5" wheels	500 lbs. (226.8 kg)	PVC	33.5"-43.25"	20.5"	1 ea
MDS86410XWWB	Steel, 5" wheels	600 lbs. (272.2 kg)	PVC	31"-41"	20"	1 ea



MDS86410XWW

## Folding paddle walkers

- Ideal design for patients with manual dexterity limitations
- The paddle mechanism can be flipped-up or pushed-down to fold the walker

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	Pkg.
MDS86410P	Aluminum, no wheels	300 lbs. (136.1 kg)	PVC	32"-39"	16.5"	4 cs
MDS86410P5	Aluminum, 5" wheels	300 lbs. (136.1 kg)	PVC	32.5"-39.5"	16.5"	4 cs
MDS86410KDBW	Aluminum, 5" wheels	350 lbs. (158.8 kg)	PVC	32.5"-39.5"	16.75"	1 ea



MDS86410P



MDS86410P5



MDS86410KDBW

\*All MDS walkers come unassembled

## Folding trigger walker

- Trigger mechanism can be pulled up to fold the walker
- Wheeled glidebrakes to help user brake easily
- Product folds for easy storage when traveling
- Comes with bag

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	Pkg.
MDS86410TRIG	Aluminum, 5" wheels	300 lbs. (136.1 kg)	PVC	29.5"-37"	17"	1 ea



MDS86410TRIG

## Guardian walkers

- All Guardian walkers come assembled
- Easy-fold sides for easy movement through narrow spaces
- Comfortable grip and no marking tips

Item No.	Description	Wt. capacity	Handle material	Height range	Hand grip width	Pkg.
G30757W	Aluminum, 5" wheels, adult	300 lbs. (136.1 kg)	PVC	32.5"-39.5"	16.5"	4 cs
G30758W	Aluminum, 5" wheels, youth	300 lbs. (136.1 kg)	PVC	27.25"-34.25"	16.5"	4 cs
G30754B	Steel, no wheels, adult	650 lbs. (294.8 kg)	PVC	32"-39"	22"	2 cs
G30755P	Aluminum, no wheels, adult	300 lbs. (136.1 kg)	PVC	32.5"-39.5"	16.5"	4 cs
G30756P	Aluminum, no wheels, youth	300 lbs. (136.1 kg)	PVC	27"-34"	16.5"	4 cs
G07749	Aluminum, no wheels, youth	200 lbs. (90.7 kg)	PVC	24.5"-28.5"	12.75"	1 ea
G07767	Aluminum, no wheels, adult	400 lbs. (181.4 kg)	PVC	32.5"-39.25"	19"	1 ea



G30757W



G30754B



G30755P



G07749

## Walker accessories

Item No.	Description	Pkg.
G07729M	7/8" footpiece, 3" wheels and glide caps	1 pr
G07731-8	1" footpiece, 3" wheels and glide caps	1 pr
G07731G	Footpiece, 3" wheels	3 pr/cs
G07722-8	Guardian, 5" front walker wheels	1 pr
MDS86615W	3" casters for walkers	2 pr/cs
MDS86615W5	Footpiece, 5" wheels and glide caps	1 pr
MDS86615W5W	Footpiece, 5" wheels and glide caps	6 pr/cs
G07724	Pediatric extension legs, 5" wheels	2 pr/cs
G07725	Walker footpiece extension kit, 5" wheels	1 kit
G07825-8	Swivel footpiece, 3" wheels	1 pr
G07722-8B	Fixed 5" bariatric dual wheels	1 ea
G07722G	Fixed 5" bariatric dual wheels	3 pr/cs
G07723M	Footpiece, 5" wheels	1 pr
MDS86615W	Footpiece, 3" wheels	1 pr
G07820-8M	For all Guardian walkers	1 pr
MDS86615SW5	Footpiece, 5" wheels	2 pr/cs
G07732G	Glidebrakes	3 pr/cs
MDS86615T	Rubber tips	1 pr
G07900CS	Glide caps	24 cs
G07904	Glide caps	6 cs
MDS86615C	1" glide caps for walkers	1 pr
MDS86615CW	1" glide caps for walkers	6 cs
MDS86615SKI	Ski gliders for walkers	6 cs
G07712	4-piece footpiece extension	1 ea
MDS86615L	10" leg replacement	4 pr/cs
G07850M	Flip tray for Guardian walkers	3 cs
G07850MB	Bariatric walker flip tray	3 cs
G07702	Walker platform attachment	1 ea
G07703M	Walker platform attachment youth	1 ea
MDS86615P	Platform attachment	1 ea
G07741	Walker front carry pouch	3 cs
G07715	Carrying basket	2 cs
G07872	Walker seat	2 cs



# Knee walkers

## 4-wheel knee walker

- Easy to maneuver with dual hand brakes for superb stopping control
- Height adjustable for patients recovering from a lower leg injury
- Can be used on right or left leg as crutch alternative

Item No.	Description	Wt. capacity	Knee pad height	Pkg.
MDS86000G4	8" wheels, Black	300 lbs. (136.1 kg)	17"-23"	1 ea



MDS86000G4

## Weil knee walker

- Tool-free, height-adjustable knee walker is a great alternative to crutches
- Universal design allows use with either a left or right below-the-knee injury for ease of use
- Features a full-length, extra-padded leg rest for maximum comfort

Item No.	Description	Wt. capacity	Knee pad height	Pkg.
MDS86000G2	8" back, 6" front wheels, Burgundy	300 lbs. (136.1 kg)	13.5"-23"	1 ea



MDS86000G2



Folds down  
for easy  
portability

# Home care beds & accessories

## Basic beds

- Headspring supports recessed to resist pinching and injuries
- Foot assembly recessed to reduce chance of pinching
- 2-year warranty on electronics, limited lifetime warranty on the frame and welds

Item No.	Description	Wt. capacity	Height range	Pkg.
MDR107002E	Semi-electric basic bed	450 lbs. (204.1 kg)	15"-20"	1 ea
MDR107003E	Full-electric basic bed	450 lbs. (204.1 kg)	15"-20"	1 ea
MDR107003ELO	Full-electric basic low bed	450 lbs. (204.1 kg)	9.5"-22.5"	1 ea



## Medlite beds

- Premium lightweight home care beds
- 5-year warranty on electronics, limited lifetime warranty on the frame and welds (all Medlite beds)

Item No.	Description	Wt. capacity	Height range	Pkg.
MDR107002L	Medlite semi-electric bed	450 lbs. (204.1 kg)	14.5"-22.5"	1 ea
MDR107003L	Medlite full-electric bed	450 lbs. (204.1 kg)	14.5"-22.5"	1 ea
MDR107003LO	Medlite full-electric low bed	450 lbs. (204.1 kg)	9.5"-22.5"	1 ea



## Bariatric bed

- Powerful motors and sturdy steel frame for heavy-duty homecare bed needs
- 2-year warranty on electronics, limited lifetime warranty on the frame and welds

Item No.	Description	Wt. capacity	Height range	Pkg.
MDR107004	Full-electric bariatric home care bed	600 lbs. (272.2 kg)	18"-26.5"	1 ea



### Side rails


- Half- and full-length Medline home care bed side rails
- 1-year warranty

Item No.	Description	Pkg.
MDR89697	Clamp on half rail	1 pr
MDR89698N	Spring loaded half rail	1 pr
MDR89694N	Spring loaded full rail	1 pr
MDS89697BAR	Bariatric half rail	1 pr
MDS89695N	Economy full rail	1 pr

MDR89697



MDR89694N




### Bed accessories

- Low bed conversion kits and transport aids for Medline home care beds

Item No.	Description	Pkg.
MDSBEDCART	Bed cart kit for home care beds	1 ea
MDSLOWBASIC	Low bed conversion kit for MDR107003E	1 ea
MDSLOWBED	Low bed conversion kit for MDR107003L	1 ea



MDSLOWBED



MDSBEDCART

### Bed assist bars

- Provides assistance getting in and out of bed with multiple hand positioning options to meet users' preferences

Item No.	Description	Handle material	Pkg.
MDS6800BA	Black	Foam	3 cs
MDS6800BAM	M-shape handle, White	Foam	1 ea
MDS6800SWIV	SwivAssist bar, Black	Foam	3 cs



MDS6800BA



MDS6800BAM



MDS6800SWIV



## Mattresses

- Cost-effective homecare bed mattresses available in three different materials with durable vinyl covers

Item No.	Description	Wt. capacity	Pkg.
MDR237827	Innerspring, 36" x 80" x 6"	250 lbs. (113.4 kg)	1 ea
MDR237827E	Fiber, 36" x 80" x 6"	250 lbs. (113.4 kg)	1 ea
MDR230981	Foam, 36" x 80" x 6"	350 lbs. (158.8 kg)	1 ea
MDR23784	Innerspring 36" x 84" x 6"	250 lbs. (113.4 kg)	1 ea
MDR23784E	Fiber, 36" x 84" x 6"	250 lbs. (113.4 kg)	1 ea
MDR211892	Foam, 42" x 84" x 6"	600 lbs. (272.2 kg)	1 ea

MDR237827



MDR237827E



MDR230981



## Overbed tables

- Standard overbed tables with various features and material options

Item No.	Description	Wt. capacity	Height range	Pkg.
MDS104015	Economy overbed table	50 lbs. (22.7 kg)	28"-45"	1 ea
MDS107015	Composite top overbed table	50 lbs. (22.7 kg)	30"-41"	1 ea
MDS108015	Fold-flat overbed table	50 lbs. (22.7 kg)	29.5"-41"	1 ea
MDS104950	Tilt-top overbed table	50 lbs. (22.7 kg)	30"-42"	1 ea

MDS104015



## Trapezes

- Free-standing and bed mounted options assist with patient transfers, movement and positioning

Item No.	Description	Wt. capacity	Pkg.
MDS80615T	Standard trapeze bar	250 lbs. (113.4 kg)	1 ea
MDS80616T	Base for MDS80615T	250 lbs. (113.4 kg)	1 ea
MDS500TPZ	Lightweight bariatric trapeze	500 lbs. (226.8 kg)	1 ea
MDS1000TPZ	Bariatric trapeze	1000 lbs. (453.6 kg)	1 ea
MDS80618T	Economy trapeze bar	250 lbs. (113.4 kg)	1 ea
MDS80618T	Economy trapeze base	250 lbs. (113.4 kg)	1 ea

MDS1000TPZ



# Canes

## Offset canes

- Single-point canes with foam handles and wrist straps

Item No.	Description	Wt. capacity	Handle material	Height range	Pkg.
MDS86420	Foam handle, Black	300 lbs. (136.1 kg)	Foam	29"-38"	6 cs
MDS86420BRZW	Foam handle, Bronze	300 lbs. (136.1 kg)	Foam	29"-38"	6 cs
MDS86420CHRW	Foam handle, Chrome	300 lbs. (136.1 kg)	Foam	29"-38"	6 cs
MDS86420FLEX	Offset smooth walker cane, Black	250 lbs. (113.4 kg)	Foam	32"-39"	1 cs



## Quad canes

- Provides a secure grip with four feet that allows greater weight-bearing and stability than single-point canes

Item No.	Description	Wt. capacity	Handle material	Height range	Pkg.
MDS86222W	Small base, Black	300 lbs. (136.1 kg)	Foam	29"-38"	2 cs
MDS86222CHR	Small base, Chrome	300 lbs. (136.1 kg)	Foam	29"-38"	2 cs
MDS86228	Large base, Black	300 lbs. (136.1 kg)	Foam	29"-38"	2 cs
MDS86228CHR	Large base, Chrome	300 lbs. (136.1 kg)	Foam	29"-38"	2 cs



## Army canes

- Medline-exclusive Army themed offset canes

Item No.	Description	Wt. capacity	Handle material	Height range	Pkg.
MDS86420AB	Black w/Army logo	300 lbs. (136.1 kg)	Foam	29"-38"	6 cs
MDS86420AC	Camo w/Army logo	300 lbs. (136.1 kg)	Foam	29"-38"	6 cs
MDS86420AG	Green w/Army logo	300 lbs. (136.1 kg)	Foam	29"-38"	6 cs



Strong in mind, body, and soul. The U.S. Army supports our nation and our values every day. Now, you can support Army service members past and present with U.S. Army licensed gear from Medline. Learn more about what it means to be a U.S. Army soldier at [www.goarmy.com](http://www.goarmy.com).



MDS86420AB



MDS86420AC



MDS86420AG

### Folding cane

- Folding capabilities make this cane perfect for daily travel

Item No.	Description	Wt. capacity	Handle material	Height range	Pkg.
MDS86430W	Black	250 lbs. (113.4 kg)	Vinyl	32"-36"	6 cs



### Standard aluminum canes

- Standard handle lets users hook the cane over their arm while opening doors or performing other tasks

Item No.	Description	Wt. capacity	Handle material	Height range	Pkg.
MDS86415	Black	300 lbs. (136.1 kg)	Foam	30"-38"	6 cs
MDS86415CHR	Chrome	300 lbs. (136.1 kg)	Foam	30"-38"	6 cs

### Bariatric canes

- Designed to give additional support that allows up to 500 lbs. weight capacity

Item No.	Description	Wt. capacity	Handle material	Height range	Pkg.
MDS86420XW	Bariatric offset, Black	500 lbs. (226.8 kg)	Foam	29"-38"	1 ea
MDS86222XW	Bariatric small base, Black	500 lbs. (226.8 kg)	Foam	29"-38"	1 ea
MDS86228XW	Bariatric large base, Black	500 lbs. (226.8 kg)	Foam	29"-38"	1 ea



MDS86415CHR



MDS86228XW

# Our complete selection will help you every step of the way.

**Medline is your trusted DME partner**, with a wide range of premium products that enhance comfort and mobility, so patients can live their best life.



**We make  
healthcare  
run better™**

**Medline Industries, Inc.**

Three Lakes Drive, Northfield, IL 60093  
Medline United States | 1-800-MEDLINE (633-5463)  
medline.com | info@medline.com

Follow us    

**Medline Canada**

1-800-268-2848 | medline.ca | canada@medline.com

**Medline México**

01-800-831-0898 | medlinemexico.com | mexico@medline.com

By federal law, licensing fees paid to the U.S. Army for the use of its trademarks provide support for the Army Trademark Licensing Program, and net licensing revenue is devoted to U.S. Army Morale, Welfare and Recreation programs. U.S. Army name, trademarks and logos are protected under federal law and used under license by [Medline Industries, Inc]. © 2021 Medline Industries, Inc. All rights reserved. Medline is a registered trademark of Medline Industries, Inc.

MKT19W841603 / LIT474 / 3M / SG / 80