


# FORCEMECH



## CARBON F1



**33** LBS\*  
POWERCHAIR WEIGHT

**300** LBS  
MAX WEIGHT CAPACITY

ELECTRIC MOTORS

**2x200**  
WATT BRUSHLESS

FULL CARBON FIBER FRAME



*Cardinal Metallic Red*

*Ocean Metallic Blue*



EASY  
FOLDING



1-877-90-FORCE

WWW.FORCEMECH.COM

\*weight without battery installed

TOTAL FREEDOM IS CALLING